
Battery Installation

- 1. Use 4 AAA batteries (included. Remove any plastic wrap from the batteries before proceeding). Always replace all 4 batteries at the same time; do not combine old and new batteries.
- 2. Open the battery compartment on the bottom of the scale.
- 3. Place the new batteries into the battery compartment according to the “+” and “-” polarity symbols at the bottom of the compartment.
- 5. Replace the battery cover.


Trouble Shooting

In Case of Difficulty...

- 1. Place your feet in the same position every time you weigh yourself. Stand still during measurement.
- 2. Weigh yourself without clothing whenever possible. Clothes can vary in weight.
- 3. Move the scale to a location where it won't be bumped. If your scale is bumped while being activated or while in use, an inaccurate reading may occur.
- 4. For the most consistent readings, weigh yourself at the same time of day. Allow for weight fluctuations due to food or drink in your system.
- 5. The maximum capacity of this scale is 440 lbs / 200 kgs / 31 st 6 lb. Should “o_Ld ” appear on the display during a weighing, this indicates that the maximum weighing capacity of the scale has been exceeded. Remove the weight immediately; otherwise, permanent damage may occur.

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
Reaching Target Weight

Maintaining proper weight is essential to healthy living. Obesity has been linked to a higher risk of heart disease, diabetes and some forms of cancer. The Body Mass Index (BMI) is the most widely accepted measurement of weight according to health professionals. Use the chart below to determine your own BMI by matching your height in the left hand column with your weight in the center. Anyone, male or female, with a BMI of 25 to 29 is considered over weight and a body mass index of 30 is considered obese. Remember to eat right and exercise.

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